



LONGER TIP
Provides a more gradual
dilation of the vessel.

ZERO-CLEARANCE TRANSITION
Between sheath and dilator, allowing
a smoother insertion into the vessel.

REGULAR AND STIFF VERSION
A steel cannula inside our stiff version provides
extra rigidity for accessing tough vessels.



MICRO-STICK® BASIC SET		
MR190421	4F MICRO-STICK®, WITH NITINOL WIRE AND 7CM ECHO TIP	10/BOX
MR190521	5F MICRO-STICK®, WITH NITINOL WIRE AND 7CM ECHO TIP	10/BOX
MR193421	4F MICRO-STICK®, WITH NITINOL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192421	4F MICRO-STICK®, WITH NITINOL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR193521	5F MICRO-STICK®, WITH NITINOL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192521	5F MICRO-STICK®, WITH NITINOL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR190422	4F MICRO-STICK®, WITH STEEL WIRE AND 7CM ECHO TIP	10/BOX
MR190522	5F MICRO-STICK®, WITH STEEL WIRE AND 7CM ECHO TIP	10/BOX
MR193422	4F MICRO-STICK®, WITH STEEL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192422	4F MICRO-STICK®, WITH STEEL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR193522	5F MICRO-STICK®, WITH STEEL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192522	5F MICRO-STICK®, WITH STEEL WIRE AND 7CM NON-ECHO TIP	10/BOX

MICRO-STICK® STIFF BASIC SET		
MR190411	4F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 7CM ECHO TIP	10/BOX
MR190511	5F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 7CM ECHO TIP	10/BOX
MR193411	4F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192411	4F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR193511	5F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192511	5F MICRO-STICK® STIFF, WITH NITINOL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR190412	4F MICRO-STICK® STIFF, WITH STEEL WIRE AND 7CM ECHO TIP	10/BOX
MR190512	5F MICRO-STICK® STIFF, WITH STEEL WIRE AND 7CM ECHO TIP	10/BOX
MR193412	4F MICRO-STICK® STIFF, WITH STEEL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192412	4F MICRO-STICK® STIFF, WITH STEEL WIRE AND 7CM NON-ECHO TIP	10/BOX
MR193512	5F MICRO-STICK® STIFF, WITH STEEL WIRE AND 4CM NON-ECHO TIP	10/BOX
MR192512	5F MICRO-STICK® STIFF, WITH STEEL WIRE AND 7CM NON-ECHO TIP	10/BOX