

## DIGNITY™ WITH MICRO-STICK® KIT

| DIGNITY™ MID-SIZED        | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------|---------------------|-----------------------|
| 9.6F Silicone Catheter    | MICTI96001          | MICTI96801            |
| 8.0F ChronoFlex® Catheter | MICTI80001          | MICTI80801            |
| 6.6F ChronoFlex® Catheter | MICTI66001          | MICTI66801            |
|                           |                     |                       |
| DIGNITY™ LOW PROFILE      | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
| 8.0F ChronoFlex® Catheter | MICTI80041          | MICTI80841            |
| 6.6F ChronoFlex® Catheter | MICTI66041          | MICTI66841            |

## DIGNITY™ WITH SILICONE FILLED SUTURE HOLES


| DIGNITY™ MID-SIZED        | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------|---------------------|-----------------------|
| 9.6F Silicone Catheter    | MRCTI9600S          | MRCTI9680S            |
| 8.0F ChronoFlex® Catheter | MRCTI8000S          | MRCTI8080S            |
| 6.6F ChronoFlex® Catheter | MRCTI6600S          | MRCTI6680S            |
|                           |                     |                       |
| DIGNITY™ LOW PROFILE      | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
| 8.0F ChronoFlex® Catheter | MRCTI8004S          | MRCTI8084S            |
| 6.6F ChronoFlex® Catheter | MRCTI6604S          | MRCTI6684S            |

## DIGNITY™ WITH SILICONE FILLED SUTURE HOLES AND MICRO-STICK® KIT

| DIGNITY™ MID-SIZED        | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------|---------------------|-----------------------|
| 9.6F Silicone Catheter    | MICTI9600S          | MICTI9680S            |
| 8.0F ChronoFlex® Catheter | MICTI8000S          | MICTI8080S            |
| 6.6F ChronoFlex® Catheter | MICTI6600S          | MICTI6680S            |
|                           |                     |                       |
| DIGNITY™ LOW PROFILE      | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
| 8.0F ChronoFlex® Catheter | MICTI8004S          | MICTI8084S            |
| 6.6F ChronoFlex® Catheter | MICTI6604S          | MICTI6684S            |

### SETS CONTAIN:

(1) Catheter (1) CT Implantable Port (2) Catheter Locks (1) Scalpel (1) 18 Ga. Introducer Needle (1) Vein Pick (1) .038" x 70cm "J" Marked Guidewire with Advancer (1) 10cc Syringe (1) Valved Tearaway Introducer (1) Tunneler (2) 22 Ga. Huber Needles (1-Straight, 1-Right Angle) (1) Blunt Tip Needle (1) Valved Tearaway Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker



### SETS WITH MICRO-STICK® CONTAIN:

(1) Catheter (1) CT Implantable Port (2) Catheter Locks (1) Scalpel (1) Vein Pick (1) .038" x 70cm "J" Marked Guidewire with Advancer (1) 10cc Syringe (1) Valved Tearaway Introducer (1) Tunneler (2) 22 Ga. Huber Needles (1-Straight, 1-Right Angle) (1) Blunt Tip Needle (1) Valved Tearaway Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker

#### ► MICRO-STICK®:

(1) 5FR Coaxial Dilator (1) .018" x 45cm Guidewire with Radiopaque Tip (1) 21 Ga. Introducer Needle with Echogenic Tip

"CHRONOFLEX® AL85AB20 (20% Barium) – main ingredient in the catheter lumen"  
"CHRONOFLEX® is a registered trademark of AdvanSource Biomaterials, Inc."



**Medical Components, Inc.**  
1499 Delp Dr., Harleysville, PA 19438, U.S.A.  
P: 215.256.4201 F: 215.256.1787  
www.medcompnet.com



MPS Medical Product Service GmbH  
Borngasse 20  
35619 Braunfels  
Germany

**ISO13485**  
Certified Quality System

CE 2797

PN2492 Rev. D 8/19

# ProFUSE<sup>CT</sup>

CT Implantable Port

## PRO-FUSE<sup>®</sup> CT LOW PROFILE WITH MICRO-STICK<sup>®</sup> KIT

| PRO-FUSE <sup>®</sup> CT LOW PROFILE  | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------------------|---------------------|-----------------------|
| 8.0F ChronoFlex <sup>®</sup> Catheter | MICTT80041          | MICTT80841            |
| 6.6F ChronoFlex <sup>®</sup> Catheter | MICTT66041          | -                     |

## PRO-FUSE<sup>®</sup> CT WITH SILICONE FILLED SUTURE HOLES


| PRO-FUSE <sup>®</sup> CT              | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------------------|---------------------|-----------------------|
| 9.6F Silicone Catheter                | MRCTT96001          | MRCTT96801            |
| 8.0F ChronoFlex <sup>®</sup> Catheter | MRCTT80001          | MRCTT80801            |
| PRO-FUSE <sup>®</sup> CT LOW PROFILE  | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
| 8.0F ChronoFlex <sup>®</sup> Catheter | MRCTT8004S          | MRCTT8084S            |
| 6.6F ChronoFlex <sup>®</sup> Catheter | MRCTT6604S          | MRCTT6684S            |

## PRO-FUSE<sup>®</sup> CT WITH SILICONE FILLED SUTURE HOLES AND MICRO-STICK<sup>®</sup> KIT

| PRO-FUSE <sup>®</sup> CT              | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
|---------------------------------------|---------------------|-----------------------|
| 9.6F Silicone Catheter                | -                   | MICTT96801            |
| 8.0F ChronoFlex <sup>®</sup> Catheter | MICTT80001          | MICTT80801            |
| PRO-FUSE <sup>®</sup> CT LOW PROFILE  | ATTACHABLE CATHETER | PRE-ATTACHED CATHETER |
| 8.0F ChronoFlex <sup>®</sup> Catheter | MICTT8004S          | MICTT8084S            |
| 6.6F ChronoFlex <sup>®</sup> Catheter | MICTT6604S          | MICTT6684S            |

### SETS CONTAIN:

- (1) Catheter (1) CT Implantable Port (2) Catheter Locks (1) Scalpel (1) 18 Ga. Introducer Needle (1) Vein Pick (1) .038" x 70cm "J" Marked Guidewire with Advancer (1) 10cc Syringe (1) Valved Tearaway Introducer (1) Tunneler (2) 22 Ga. Huber Needles (1-Straight, 1-Right Angle) (1) Blunt Tip Needle (1) Valved Tearaway Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker



### SETS WITH MICRO-STICK<sup>®</sup> CONTAIN:

- (1) Catheter (1) CT Implantable Port (2) Catheter Locks (1) Scalpel (1) Vein Pick (1) .038" x 70cm "J" Marked Guidewire with Advancer (1) 10cc Syringe (1) Valved Tearaway Introducer (1) Tunneler (2) 22 Ga. Huber Needles (1-Straight, 1-Right Angle) (1) Blunt Tip Needle (1) Valved Tearaway Introducer IFU (1) Port IFU (1) Patient Information Pack (1) Patient Chart Sticker

#### ► MICRO-STICK<sup>®</sup>:

- (1) 5FR Coaxial Dilator (1) .018" x 45cm Guidewire with Radiopaque Tip (1) 21 Ga. Introducer Needle with Echogenic Tip



"CHRONOFLEX<sup>®</sup> AL85AB20 (20% Barium) – main ingredient in the catheter lumen"  
 "CHRONOFLEX<sup>®</sup> is a registered trademark of AdvanSource Biomaterials, Inc."