



#### HEMO-CATH® LT

**Indications for Use:** The Hemo-Cath® Dialysis Catheter is indicated for use in attaining Long-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required. The curved Hemo-Cath® Catheter is intended for internal jugular vein insertion.

**Contraindications:** This catheter is intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions.


Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.




<b>Silicone</b> Material	<b>Step Tip</b> Tip Design	<b>8F, 10F, 12.5F</b> French Sizes	<b>Straight, Pre-Curved</b> Configurations	<b>Antegrade</b> Insertion Type
-----------------------------	-------------------------------	---------------------------------------	---	------------------------------------



CIRCLE "C" INTERNAL LUMEN DESIGN



LONGER LENGTHS



SILICONE MATERIAL



ANGLED ARTERIAL TIP



STRAIGHT CATHETER BASIC SET		
SL18P	8F X 18CM SILICONE HEMO-CATH® BASIC SET	5/BOX
SL24P	8F X 24CM SILICONE HEMO-CATH® BASIC SET	5/BOX

STRAIGHT CATHETER BASIC SET		
MC101041P	10F X 15CM SILICONE HEMO-CATH® LT	5/BOX
MC101042P	10F X 18CM SILICONE HEMO-CATH® LT	5/BOX
MC101043P	10F X 20CM SILICONE HEMO-CATH® LT	5/BOX
MC101044P	10F X 22CM SILICONE HEMO-CATH® LT	5/BOX
MC101045P	10F X 24CM SILICONE HEMO-CATH® LT	5/BOX

STRAIGHT CATHETER BASIC SET		
MC101241	12.5F X 15CM SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX
MC101242	12.5F X 18CM SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX
MC101243	12.5F X 24CM SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX

STRAIGHT CATHETER BASIC SET		
SL28E.	12.5F X 28CM SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX
SL32E.	12.5F X 32CM SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX

PRE-CURVED CATHETER BASIC SET		
SL28PCE.	12.5F X 28CM PC SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX
SL32PCE.	12.5F X 32CM PC SILICONE HEMO-CATH® DOUBLE LUMEN	5/BOX

**8F BASIC SET CONTAINS:**

- (1) Catheter
- (1) Tunneler
- (1) Peelable Introducer
- (1) Guidewire
- (1) 10cc Syringe
- (1) Adhesive Wound Dressing
- (1) Catheter Clip
- (1) Introducer Needle
- (2) End Caps

**10F SET CONTAINS:**

- (1) Catheter
- (1) Introducer Needle
- (1) Adhesive Wound Dressing
- (1) J/Flex Guidewire
- (1) Tunneler
- (1) Catheter Clip
- (2) End Caps
- (1) 10cc Syringe
- (1) Peelable Introducer

**12.5F x 15, 18, 24CM SET CONTAINS:**

- (1) Catheter
- (2) End Caps
- (1) Scalpel
- (1) Peelable Introducer
- (1) Introducer Needle
- (1) Guidewire
- (1) Adhesive Wound Dressing
- (1) Catheter Clip
- (1) Tunneler

**12.5F x 28, 32CM SET CONTAINS:**

- (1) Catheter
- (2) End Caps
- (1) Peelable Introducer
- (1) Introducer Needle
- (1) Guidewire
- (1) Tunneler
- (1) Catheter Clip
- (Straight Catheter Only)

Refer to the Table of Contents for Important Risk Information regarding this device.