



SUBCLAVIAN

Indications for Use: The Subclavian Catheter is designed for acute hemodialysis and apheresis. It may be inserted percutaneously and the preferred site is the subclavian vein. This catheter is indicated for a duration less than (30) days.

Contraindications: This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

SUBCLAVIAN

Polyurethane

Material

Tapered Tip


Tip Design

8F

French Size

**Single Lumen,
Straight**

Configurations



ROTATING SUTURE WING

REMOVABLE "Y" ADAPTOR
And pre-loaded stylet.

CATHETER ONLY		
MCY306PS	8F X 15CM SUBCLAVIAN CATHETER, WITH REMOVABLE "Y" HUB	10/BOX
MCY308PS	8F X 20CM SUBCLAVIAN CATHETER, WITH REMOVABLE "Y" HUB	10/BOX

BASIC SET		
MCYK306PSE*	8F X 15CM SUBCLAVIAN CATHETER, WITH REMOVABLE "Y" HUB	10/BOX
MCYK308PSE*	8F X 20CM SUBCLAVIAN CATHETER, WITH REMOVABLE "Y" HUB	10/BOX
MCB304K	8F X 12CM SINGLE LUMEN CATHETER KIT, PRE-CURVED	10/BOX
MCB306K	8F X 15CM SINGLE LUMEN IJ CATHETER KIT, PRE-CURVED	10/BOX
MCB308K	8F X 20CM SINGLE LUMEN IJ CATHETER KIT, PRE-CURVED	10/BOX

CATHETER ONLY SET CONTENTS:

- (1) Catheter w/Stylet
- (1) Removable "Y" Adaptor
- (1) Injection Cap
- (1) Patient I.D. Card

BASIC SET CONTAINS:

- (1) Catheter w/ Stylet
- (1) Removable "Y" Adaptor*
- (2) End Caps
- (1) Introducer Needle
- (1) Guidewire
- (1) Cap
- (1) Adhesive Wound Dressing

[Refer to the Table of Contents for Important Risk Information regarding this device.](#)