



T-3[®] CT

Indications for Use: The Medcomp[®] T-3[®] CT Catheter is a triple lumen catheter indicated for use in attaining short-term vascular access for hemodialysis, apheresis. The third internal lumen is intended for infusion, power injection of contrast media and central venous pressure monitoring. The catheter is intended to be inserted in the jugular, femoral or subclavian vein as required. The maximum recommended infusion rate is 5ml/sec for power injection of contrast media.


Contraindications: This catheter is intended for short-term (less than 30 days) vascular access only and should not be used for any purpose other than indicated in these instructions. This device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient to accommodate the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors may prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



TRIPLE INTERNAL LUMEN DESIGN

Polyurethane Material	Step Tip Tip Design	15.5F French Size	Straight Configuration
---------------------------------	-------------------------------	-----------------------------	----------------------------------



CT CATHETER ONLY		
MC013541	15.5F X 15CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013542	15.5F X 20CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013543	15.5F X 24CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013544	15.5F X 28CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013545	15.5F X 32CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX

CT CATHETER ONLY CONTAINS:
 (1) Catheter
 (3) Dilators
 (1) Valved Peelable Introducer
 (1) Patient ID Card
 (3) End Caps

CT BASIC SET		
MC013501	15.5F X 15CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013502	15.5F X 20CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013503	15.5F X 24CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013504	15.5F X 28CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX
MC013505	15.5F X 32CM T-3 [®] CT, STRAIGHT EXTENSION	5/BOX

CT BASIC SET CONTAINS:
 (1) Catheter
 (3) Dilators
 (1) Scalpel
 (1) Valved Peelable Introducer
 (1) Guidewire
 (1) Adhesive Wound Dressing
 (1) 2-0 Silk Suture
 (1) Patient ID Card
 (3) End Caps
 (1) Introducer Needle

Refer to the Table of Contents for Important Risk Information regarding this device.